



东莞市弗瑞斯生物科技有限公司
Dongguan Fores Biotechnology Co., Ltd.

4th Floor, No.2 Building, No.5, Shipai Industrial Avenue, Shipai Town, Dongguan City,
Guangdong Province, China

Tel: +86-0769-82203796

Email: liang@choebe.com

Product name	Disposable Face Mask
Product number	Y001
Sizes	175*95mm
Photo	
Application	Used in laboratory, food processing, beauty salon, electronics industry, etc.
Material	<ol style="list-style-type: none"> 1. Spun-bonded non-woven fabric 2. Melt-blown non-woven fabric 3. ear loops and nose bridge bar 4. 3-layer Fiber Gabric Construction
Fabric Density	<p>3 layer:</p> <p>Inner layer : 20 g/m2 PP non-woven fabrics,</p> <p>Middle filter : 25 g/m2 melt blown filter</p> <p>Outer layer : 25g/m2 PP non-woven fabrics</p>
Features	<ol style="list-style-type: none"> 1. Three-dimensional nose bridge strip for easy fit 2. Extra-soft ear loops and flexible binding bridge 3. Multilayered Protection and Extra breathability
Instruction	<ol style="list-style-type: none"> 1. Flatten the mask, wear white color inside, put nosepiece in the top area and place loops over ears. 2. Adjust nosepiece for dust protection. 3. Pull pleats down below chin to get breathing space.
Packing	50pcs in a polybag, 1 bags/inner box, 40boxes/carton, 2000pcs/carton, G. W. :8.2 KGS, carton size:515*395*340mm



东莞市弗瑞斯生物科技有限公司
Dongguan Fores Biotechnology Co., Ltd.

4th Floor, No.2 Building, No.5, Shipai Industrial Avenue, Shipai Town, Dongguan City,
Guangdong Province, China

Single use or not	Single use
Sterile or not	Non sterilized
Attention	Once the face mask gets dirty and can not provide further protection, please change another new one
Storage	Store in a dry, cool place (<40° C) away from direct sunlight and fluorescent light