

东莞市弗瑞斯生物科技有限公司

Dongguan Fores Biotechnology Co., Ltd.

No.5, Shipai Industrial Avenue, Shipai Town, Dongguan City, Guangdong Province, China

Tel: +86-0769-82203796 Email: <u>liang@choebe.com</u>

Product name	KN95 Protective Mask
Product number	K001
Sizes	170*100mm
Photo	RN95
Application	Mask for protection use only
Material	 Spun-bonded non-woven fabric Melt-blown non-woven fabric ear loops and nose bridge bar 5-layer Construction
Fabric Density	5 layer: First layer: 50g non-woven fabrics, Second layer: 50g heat sealing cotton Third layer: 25g melt blown filter Fourth layer: 25g melt blown filter Fifth layer: 30g non-woven fabrics
Instruction &Cautions	 Please do not use and promptly replace the mask if the package and mask are damaged or stained. Please read the instructions and directions carefully before use, improper usage may cause the product's protection performance to fail. Prohibited to use after expiration date. It is not recommended to use the same mask for a long time in non-special circumstances. Please replace mask at least in a week.



东莞市弗瑞斯生物科技有限公司

Dongguan Fores Biotechnology Co., Ltd.

4thr Floor, No. 2 Building, No. 5, Shipai Industrial Avenue, Shipai Town, Dongguan City, Guangdong Province, China

Packing	1pcs in a polybag, 50 bags/inner box, 20boxes/carton, 1000pcs/carton, G.W.:9.60 KGS,
	carton size:59.5*32.5*42.5cm
Single use or	Single use
not	
Sterile or not	Non sterilized
Attention	Once the face mask gets dirty and can not provide further protection, please change another new one
Storage	10-25°C temperature, relative humidity <80%, keep cool and ventilated